



Größe	1	2	3	4	5	6	7	8
Fußlänge A/cm	22	23	24	25	26	27	28,5	29,5
Ballenbreite/mm	75	80	85	90	95	95	100	100
Absatzhöhe B/mm	15	15	15	15	15	15	15	15
Fußhöhe C/mm	75	75	75	80	80	80	90	90
Gewicht/g	230	260	280	330	360	400	440	500

REF

21760...	links	122	123	124	125	126	127	128	129
	rechts	222	223	224	225	226	227	228	229

passend für SACH-Fuß-Adapter M10